

ActiPix D200

UV-VIS Imaging Light Box

Flexible Wavelength Options

Designed with stability and flexibility in mind the ActiPix™ D200 series Light Box can be configured to suit user needs.

Light Box Specifications

- LED wavelengths : 255 / 280 / 320 / 525 +/-5 nm*
- Pulsed xenon flash lamp wavelengths 180–800 nm via bandpass filter (214 nm default)
- LED 20 mA power
- Xenon Lamp 6 W, 0.27 μ F
- Xenon expected life time $>1 \times 10^9$ flashes
- LED expected 1000-5000 hours
- Light Box weight 400 –700 g (configuration dependent)
- Filter kit for narrow bandwidth illumination available
- Power 24V DC via H200 Control Hub

*others available



For more information contact:

Paraytec Ltd
York House, Outgang Lane,
Osbaldwick, York, YO19 5UP
info@paraytec.com
+44 (0)1904 436620