

ActiPix D200

UV-VIS Imaging Control Hub

Flexible Control Options

The ActiPix™ D200 Control Hub coordinates the sequencing of the various components of the D200 Imaging Platform and the power distribution.

Control Hub Specifications

Neatly couples to the base of the detector module or can be used remotely with optional extra length cables

- Control Hub dimensions (L x W x H)
160 x 100 x 30 mm
- Control Hub weight 360 g
- Power input 100 – 240 V AC
50/60 Hz via external power supply
- Power output 4 x 24 V DC
- USB In 2 x USB 2.0 / Out 1 x USB 2.0
- Certified to IEC 61010



For more information contact:

Paraytec Ltd
York House, Outgang Lane,
Osbaldwick, York, YO19 5UP
info@paraytec.com
+44 (0)1904 436620

